

# Quick Start Guide

This guide will help you prepare the Oculus **Quest 1** (black) or **Quest 2** (white) and get you into the PrecisionOS VR modules. To get started, you will need the following: an Oculus headset, two controllers, a Wi-Fi network, and an obstruction free area (6.5ft x 6.5ft recommended).

---

## 1. POWER UP THE HEADSET

Press the Power button for 2–3 seconds, until you see the Oculus logo.

Note the slightly different button locations between the Quest 1 and the Quest 2, respectively.



---

## 2. PICK UP BOTH CONTROLLERS

Hold in pistol-grip with your index fingers on the Trigger buttons (circled). Pull either Trigger buttons to interact with elements in VR.

Note: The Oculus Home button is on the **right** controller and is used to open the Oculus dashboard anywhere in VR – use it to access the device settings, exit an app, launch an app, take screenshots/videos, etc.



---

## 3. PUT THE HEADSET ON

Pull and position the headset until it sits comfortably, and the side straps cradle the ears. You may need to manually adjust the Velcro straps.



---

## 4. CREATE A “GUARDIAN” BOUNDARY


Follow the in-display instructions to calibrate the floor height and define a safe play area; “room scale” is recommended for a more immersive experience.

Note: Scan this QR code for an instructional video if you need help with this step.




## 5. CONNECT TO WI-FI NETWORK

Press the Oculus Home button on the controller and select the Wi-Fi icon (circled below) using the Trigger. Aim your cursor at the network list and use the thumbstick on the controller to scroll up and down; select your network and enter your password, then select “Go”.




## 6. INSTALL THE PRECISION OS APPS

From the Oculus dashboard, select the “Apps” icon, or “Library” icon if on the older dashboard (circled in red below). Then, change the dropdown filter in the corner to “Not Installed” (circled in yellow below) to see the PrecisionOS apps. Select each app with the Trigger to install.



## 7. OPEN THE PRECISION OS LAUNCHER


While still in the same Apps screen as in Step 6, change the same dropdown filter back to “All”. Using the thumbstick to scroll up and down if needed, locate the “Platform Launcher” app and select it with the trigger. Installing the “Platform Launcher” app is **required** in order for the other PrecisionOS apps to work.



---

## 8. SIGN IN, REGISTER OR GUEST MODE

Using your now virtually gloved hands to tap the virtual keys (no triggers needed), sign-in, register a free account with your email, or select the “Guest” option. You can access the modules in Guest mode but you won't be able to see your metrics on the web-based dashboard.



---

## 9. EXPLORE THE PROCEDURES

Select a procedure by “virtually” tapping it with your gloved hand. Once in the OR, you can follow the written and spoken instructions from the virtual assistant to complete the procedure. Tap “Home” and then “Sign Out” when finished.

---

## FREQUENTLY ASKED QUESTIONS

**My Oculus app is asking for a 6-digit code to confirm my login:** locate the serial number on your headset, then call PrecisionOS Support at 1-877-673-0176 (ext. 1).

**I don't see my hands or controllers in VR:** check the batteries in the controllers, replace as needed, then wipe the camera sensors on the front of the headset with a microfiber cloth.

**If you encounter any issues or otherwise need assistance with the experience, please contact PrecisionOS Support at [support@precisionostech.com](mailto:support@precisionostech.com), or call 1-877-673-0176 (ext. 1).**