




OCULUS QUEST II HANDSET CHEAT SHEET

#1 Index Finger / Front "TRIGGERS": Picking Up / Gripping Objects

Left Side Buttons
Similar to Right Side

Left Flat Button:
NOT USED (YET)



"B" Button:
Take an X-Ray ("host" only)
(Distal - Raised)

"Joystick" (Gestures):
Push UP = Index Finger Point
Pull DOWN = Thumbs up

#3 "OCULUS" Button:
Exit App / "Esc" Key
(Flat - Proximal to Joystick)

#2 "A" Button:
Teleporting
(Hold, Aim, Release)
(Proximal - Raised)

Side Gripping "Triggers":
NOT USED (YET)