

AlloAid SubChon™

ENGINEERED ALLOGRAFT



Tapered Design Resists Subsidence

Osteoconductive Porous Scaffold

Size Specific Instrumentation




AlloAid SubChon™

ENGINEERED ALLOGRAFT




ORDER INFO

| CAT NO. | DESCRIPTION | DIA (A) X LENGTH (B) |
|---------|--------------------------------|----------------------|
| ASB 005 | AlloAid SubChon Implant - 5mm | 5 x 10mm |
| ASB 007 | AlloAid SubChon Implant - 7mm | 7 x 10mm |
| ASB 009 | AlloAid SubChon Implant - 9mm | 9 x 10mm |
| ASB 011 | AlloAid SubChon Implant - 11mm | 11 x 11mm |



AlloAid SubChon Instrument Kits

| CAT NO. | DESCRIPTION |
|-----------|--|
| C01 S1001 | AlloAid SubChon Sizing Kit, Single-Use, Sterile |
| C01 S1505 | AlloAid SubChon Instrument Kit, Size 5mm, Single-Use, Sterile |
| C01 S1507 | AlloAid SubChon Instrument Kit, Size 7mm, Single-Use, Sterile |
| C01 S1509 | AlloAid SubChon Instrument Kit, Size 9mm, Single-Use, Sterile |
| C01 S1511 | AlloAid SubChon Instrument Kit, Size 11mm, Single-Use, Sterile |



AlloAid StimDisk®

This membrane consists of fibrillar and membranous collagens, elastin, and a mix of growth factors and cytokines.



ORDER INFO

| CAT NO. | DESCRIPTION | DIAMETER | CORRESPONDING SUBCHON SIZE |
|---------|------------------|----------|----------------------------|
| ASD 006 | AlloAid StimDisk | 6mm | 5mm |
| ASD 008 | AlloAid StimDisk | 8mm | 7mm |
| ASD 010 | AlloAid StimDisk | 10mm | 9mm |
| ASD 012 | AlloAid StimDisk | 12mm | 11mm |

All content contained herein is furnished for informational purposes only. In2Bones does not recommend a particular surgical product or procedure suitable for all patients. Each surgeon must evaluate the appropriateness of a device and corresponding techniques based on medical training, clinical judgment and surgical experience. The proper surgical technique and/or procedure are the responsibility of the medical professional. Indications, contraindications, warnings, and precautions are listed in the implant package insert and should be reviewed carefully by the physician and operating room personnel prior to any proposed procedure. Availability of these products might vary from a given country or region to another as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of a physician.



Corporate Headquarters

In2Bones Global, Inc. • Memphis, TN • USA
844. 602. 6637 • Info@i2b-USA.com

International Office

In2Bones SAS • Lyon • France
+33 (0)4 72 29 26 26

FOLLOW US....



In2Bones.com

™ and © denote trademarks and registered trademarks of In2Bones or its affiliates. In2Bones USA Memphis TN, 38119 USA / In2Bones SAS, 69130 Ecully, France © 2022 In2Bones USA, Memphis, TN • Patent Pending • All rights reserved • SP-21005-02 Rev. B, 1222