



## Disposable Monopolar Scissors



A universal style handle, precise cutting and sustained sharpness make these disposable scissors clinically superior and cost effective.



**LAPAROSCOPIC  
SPECIMEN BAG**



**DETACHATIP® II  
SYSTEM**



**REFLEX® ELC 530  
AUTOMATIC CLIP APPLIER**



To learn more about these and other innovative products, call **800-448-6506** or visit **ConMed.com**.



525 French Road, Utica, NY 13502

These Disposable Monopolar Scissors offer precise and smooth cutting – all the way to the tips. Distal tip sharpness provides exceptional performance for partial opening of cystic duct in preparation for cholangiogram.

### Comfortable Design

With effective one hand/finger operation and scissor blade orientation control, the ergonomically enhanced handle delivers a comfortable tactile feel and control for all hand sizes. The rotation knob allows easy 360° rotation of the shaft.

### Reliable Operation

Extended shaft insulation to the blade minimizes the chance of thermal damage through stray burns from the shaft. And with no wings on open blades, the risk of breaching insulation is also reduced.

**Precise and sharp blades provide smooth cutting especially at blade tips**

**Universal 4mm monopost fits most cords**

**Extended PTFE insulation to the blade reduces chance of thermal damage**

**Internal hinge design for open blades reduces risk of breaching insulation**



## Disposable Monopolar Scissors Ordering Information

Product	Diameter	Length	Packaging	Catalog Number
Mini Metz Disposable Scissors	5mm	33cm	10/box	1221
Metz Disposable Scissors	5mm	33cm	10/box	1222

Customer Service: 1800-448-6506  
International Sales: +1(315) 797-8375  
Fax: 1800-438-3051

Website: [www.conmed.com](http://www.conmed.com)  
Email: [info@mail.conmed.com](mailto:info@mail.conmed.com)